


Product Name: Recombinant Human DHX38 Protein
Catalog #: PEH5193



Summary

Name	PRP16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92620
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RMLTSGVVHR LEVDEDFEED NAAKVHLMVH NLVPPFLDGR IVFTKQPEPV IPVKDATSDL AIIARKGSQT VRKHREQKER KKAQHKHWEL AGTKLGDIMG VKKEEEDPKA VTEDGKVDYR TEQKFADHMK RKSEASSEFA

image



Background

Alternative Names Pre-mRNA-splicing factor ATP-dependent RNA helicase PRP16 (EC 3.6.4.13) (ATP-dependent RNA helicase DHX38) (DEAH box protein 38)

Note

For Research Use Only , Not for Diagnostic Use