


Product Name: Recombinant Human SORL1 Protein
Catalog #: PEH5194



Summary

Name	SORL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92673
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PNTIYIERHE PSGYSTVFRS TDDFFQSRENQ EVILEEVRDF QLRDKYMFAT KVVHLLGSEQ QSSVQLWVSF GRKPMRAAQF VTRHPINEYY IADASEDQVF VCVSHSNNRT NLYISEAEGF KFSLSLENVL YYSPGGAGSD TLVRYFANEP FADFHRVEGL QGVYIATLIN GSMNEENMRS VITFDKGGTW EFLQAPAFGT YGEKINCELS QGCSLHLAQR LSQLLNLQLR RMPILSKESA PGLIATGSV GKNLASKTNV YISSAGARW REALPGPHY TWGDHGGIIT AIAQGMETNE LKYSTNEGET WKTFIFSEKP VFVYGLLTP GEKSTVFTIF GSNKENVHSW LILQVNATDA LGVPCTENDY KLWSPSDEGR NECLLGHKTV FKRRTPHATC FNGEDFDRPV VVSNCSTRE DYECDFGFKM SEDLSLEVCV PDPEFSGKSY SPPVPCPVGS TYRRTGRYRK ISGDTCSGGD VEARLEGELV PCPLAEENEF

image



Background

Product Name: Recombinant Human SORL1 Protein
Catalog #: PEH5194



Alternative Names

Sortilin-related receptor (Low-density lipoprotein receptor relative with 11 ligand-binding repeats) (LDLR relative with 11 ligand-binding repeats) (LR11) (SorLA-1) (Sorting protein-related receptor containing LDLR class A repeats) (SorLA)

Note

For Research Use Only , Not for Diagnostic Use