


Product Name: Recombinant Human RSC1A1 Protein
Catalog #: PEH5195



Summary

Name	RSCA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92681
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NHPARSSGQS PDVGNPMSLA RSVSASVCPI KPSDSRIEP KAVKALKASA EFQLNSEKKE HSLQDLSDH ASSADHAPTD QSPAMPMQNS SEEITVAGNL EKSAERSTQG LKFHLHTRQE ASLSVTSTRM HEPQMFLGEK DWHPENQNLS QVSDPQQHEE PGNEQYEVAQ QKASHDQEYL CNIGDLELPE ERQQNQHKIV DLEATMKGNG LPQNVDPPSA KKSIPSECS GCSNSETFME IDTAQQSLVT

image



Background

Alternative Names	Regulatory solute carrier protein family 1 member 1 (Transporter regulator RS1) (hRS1)
--------------------------	--

Note

Product Name: Recombinant Human RSC1A1 Protein
Catalog #: PEH5195



For Research Use Only , Not for Diagnostic Use