


Product Name: Recombinant Human RREB1 Protein
Catalog #: PEH5196



Summary

Name	RREB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92766
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFKLAEGDGE AGAGGAASQE QKLACDTCGK SFKFLGTLNR HRKAHGRQEP KDEKGDGAST AEEGPQPAPE QEEKPPETPA EVVESAPGAG EAPAEKLAEE TEGPSDGESA AEKRSSEKSD DDKKPKTDSP KSVASKADKR KKVCSVCNKR FWSLQDLTRH MRSHTGERPY KCQTCERTFT LKHS�VRHQR

image



Background

Alternative Names	Ras-responsive element-binding protein 1 (RREB-1) (Finger protein in nuclear bodies) (Raf-responsive zinc finger protein LZ321) (Zinc finger motif enhancer-binding protein 1) (Zep-1)
--------------------------	--

Note

Product Name: Recombinant Human RREB1 Protein
Catalog #: PEH5196



For Research Use Only , Not for Diagnostic Use