


Product Name: Recombinant Human FOXJ1 Protein
Catalog #: PEH5198



Summary

Name	FOXJ1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92949
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALPPGGTDPH GYHQVPGSAA PGSPLAADPA CLGQPHTPGK PTSSCTSRSA PPGLQAPPPD DVDYATNPHV KPPYSYATLI CMAMQASKAT KITLSAIYKW ITDNFCYFRH ADPTWQNSIR HNLSLNKCFI KVPREKDEPG KGGFWRIDPQ YAERLLSGAF KKRRLPPVHI HPAFARQAAQ EPSAVPRAGP LTVNTEAQQ LREFEEATGE AGWGAGEGRL GHKRRKQLPK RVAKVPRPPS TLLPTPEEQG

image



Background

Alternative Names Forkhead box protein J1 (Forkhead-related protein FKHL13) (Hepatocyte nuclear factor 3 forkhead homolog 4) (HFH-4)

Note

Product Name: Recombinant Human FOXJ1 Protein
Catalog #: PEH5198



For Research Use Only , Not for Diagnostic Use