


Product Name: Recombinant Human CREB3L1 Protein
Catalog #: PEH5206



Summary

Name	CR3L1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96BA8
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LVQMPPTPPS SHGSDSDGSQ SPRSLPPSSP VRPMARSSTA ISTSPLLTA P HKLQGTSGPL LLTEEEKRTL IAEGYPIPTK LPLTKAEEKA LKRVRRIKIN KISAQESRRK KKEYVECLEK KVETFTSENN ELWKKVETLE NANRTLLQQL

image



Background

Alternative Names	Cyclic AMP-responsive element-binding protein 3-like protein 1 (cAMP-responsive element-binding protein 3-like protein 1) (Old astrocyte specifically-induced substance) (OASIS) [Cleaved into: Processed cyclic AMP-responsive element-binding protein 3-like protein 1]
--------------------------	---

Note

Product Name: Recombinant Human CREB3L1 Protein
Catalog #: PEH5206



For Research Use Only , Not for Diagnostic Use