


Product Name: Recombinant Human ZNF653 Protein
Catalog #: PEH5212



Summary

Name	ZN653_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96CK0
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGSQPSTMDA TAVAGIETKK EKEDLCLLKK EEKEEPVAPE LATTVPESAE PEAEADGEEL DGSDMSAIY EIPKEPEKRR RSKRSRVMDA DGLLEMFHCP YEGCSQVYVA LSSFQNHVNL VHRKGKTKVC

image



Background

Alternative Names Zinc finger protein 653 (67 kDa zinc finger protein) (Zinc finger protein Zip67)

Note

For Research Use Only , Not for Diagnostic Use