


Summary

Name	MTMR8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96EF0
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QSMLESLEI KKQRAMLETD VHELEKCLKV RDEPPEEICT CSQLGNILSQ HLGSPLTNPL GFMGINGDLN TLMENGTLSR EGGLRAQMDQ VKSQGADLHH NCCEIVGSLR AINISGDVGI SEAMGISGDM CTFEATGFSK DLGICGAMDI SEATGISGNL GISEARGFSG DMGILGDTGI SKASTKEADY

image



Background

Alternative Names	Myotubularin-related protein 8 (Phosphatidylinositol-3,5-bisphosphate 3-phosphatase) (EC 3.1.3.95) (Phosphatidylinositol-3-phosphate phosphatase) (EC 3.1.3.64)
--------------------------	---

Note

Product Name: Recombinant Human MTMR8 Protein
Catalog #: PEH5215



For Research Use Only , Not for Diagnostic Use