


**Product Name: Recombinant Human TCEAL4 Protein**  
**Catalog #: PEH5216**



## Summary

|                              |                                                                                                                                                                                                                                                          |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | TCAL4_HUMAN                                                                                                                                                                                                                                              |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                      |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                         |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                    |
| <b>Accession #</b>           | Q96EI5                                                                                                                                                                                                                                                   |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                   |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                    |
| <b>Molecular Mass</b>        | 25 kDa                                                                                                                                                                                                                                                   |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                       |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                           |
| <b>Sequence</b>              | MEKLYSENEG MASNQGKMEN EEQPQDERKP EVTCTLEDKK LENEGKTENK<br>GKTGDEEMLK DKGKPESEGE AKEGKSEREG ESEMEGGSER EGKPEIEGKP<br>ESEGEPGSET RAAGKRPAED DVPRKAKRKT NKGLAHYLKE YKEAIHDMNF<br>SNEDMIREFD NMAKVQDEKR KSKQKLGAFL WMQRNLQDPF YPRGPREFRG<br>GCRAPRRDIE DIPYV |

## image



## Background

|                          |                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Transcription elongation factor A protein-like 4 (TCEA-like protein 4)<br>(Transcription elongation factor S-II protein-like 4) |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------|

## Note

**Product Name: Recombinant Human TCEAL4 Protein**  
**Catalog #: PEH5216**



---

For Research Use Only , Not for Diagnostic Use