



Product Name: Recombinant Human DAZAP1 Protein
Catalog #: PEH5217



Summary

Name	DAZP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96EP5
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLFVGGLDWS TTQETLSYF SQYGEVDCV IMKDKTTNQS RGFGFVKFKD PNCVGTVLAS RPHTLDGRNI DPKPCTPRGM QPERTRPKEG WQKGPRSDNS KSNKIFVGGI PHNCGETELR EYFKKFGVVT EVVMYDAEK QRPRGFGFIT FEDEQSVDAQ VNMHFHDIMG KKVEVKRAEP

image



Background

Alternative Names DAZ-associated protein 1 (Deleted in azoospermia-associated protein 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DAZAP1 Protein
Catalog #: PEH5217

