


Product Name: Recombinant Human SIPA1 Protein
Catalog #: PEH5219



Summary

Name	SIPA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96FS4
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SFLPRTLRLR NSISIRIMSEA GSGTLEDEWQ AISEIASTCN TILESLSREG QPIPEGDPK GTPKSDAEPE PGNLSEKVSH LESMLRKLQE DLQKEKADRA

image



Background

Alternative Names	Signal-induced proliferation-associated protein 1 (Sipa-1) (GTPase-activating protein Spa-1) (p130 SPA-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use