


Product Name: Recombinant Human DDRGK1 Protein
Catalog #: PEH5224



Summary

Name	DDRGK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96HY6
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGRPRRRDL GSRLQAQRRR QRVAAEAE NEEAVILAQ EEGVEKPAE THLSGKIGAK KLRKLEEKQA RKAQREAEER EREERKRLES QREAEWKKEE ERLRLEEEQK EEEERKAREE QAQREHEEYL KLKEAFVVEE EGVGETMTEE QSQSFLTEFI NYIKQSKVVL LEDLASQVGL RTQDTINRIQ DLLAEGTITG VIDDRGKFIY ITPEELAAVA NFIRQRGRVS IAELAQASNS

image



Background

Alternative Names	DDRGK domain-containing protein 1 (Dashurin) (UFM1-binding and PCI domain-containing protein 1)
--------------------------	---

Note

Product Name: Recombinant Human DDRGK1 Protein
Catalog #: PEH5224



For Research Use Only , Not for Diagnostic Use