


Product Name: Recombinant Human ZNF625 Protein
Catalog #: PEH5225



Summary

Name	ZN625_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96I27
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KPYECKQCGK AFSHSGSLRI HERTHTGEKP YECSECGKAF HSSTCLHAHK ITHTGEKPYE CKQCGKAFVS FNSVRYHERT HTGEKPYECK QCGKAFRSAS HLRTHGRTHT GEKPYECKQC GKAFGCASSV KIHERTHTGE KPCSSNTSKG

image



Background

Alternative Names Zinc finger protein 625

Note

For Research Use Only , Not for Diagnostic Use