


Product Name: Recombinant Human ZNF594 Protein
Catalog #: PEH5228



Summary

Name	ZN594_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96JF6
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEEKKSARAA SEKLQRQITQ ECELVETSNS EDRLKHWVS PLKDAMRHLP SQESGIREMH IIPQKAIVGE IGHGCNEGEK ILSAGESHR YEVSQGNFKQ KSGLTEHQKI HNINKTYECK

image



Background

Alternative Names Zinc finger protein 594 (Zinc finger protein HZF18)

Note

For Research Use Only , Not for Diagnostic Use