


Product Name: Recombinant Human ZNF462 Protein
Catalog #: PEH5231



Summary

Name	ZN462_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96JM2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VVPIEVCRSK LSKYLQGVVF RCDKCTFTCS SDESLQQHIE KHNELKPYKC QLCYETKHT EELDSHLRDE HKVSRNFELV GRVNLQLEQ MKEKMESSSS DDEDKEEEMN SKAEDRELMR FSDHGAALNT EKRFPCFCG RAFSQGSEWE RHVLRHGMAL NDTKQVSREE IHPKEIMENS VKMPSIEEKE DDEAIGIDFS LKNETVAICV VTADKSLLEN

image



Background

Alternative Names Zinc finger protein 462 (Zinc finger PBX1-interacting protein) (ZFPIP)

Note

Product Name: Recombinant Human ZNF462 Protein
Catalog #: PEH5231



For Research Use Only , Not for Diagnostic Use