


Product Name: Recombinant Human CCDC136 Protein
Catalog #: PEH5232



Summary

Name	CC136_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96JN2
Host	E.coli
Species	Human
Molecular Mass	69 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSLTETELEE LRAQVLQLVA ELEETRELAG QHEDDSLELQ GLLIEDERLAS AQQAEVFTKQ IQQLQGELRS LREEISLLEH EKESELKEIE QELHLAQAEI QSLRQAAEDS ATEHESDIAS LQEDLCRMQN ELEDMERIRG DYEMEIASLR AEMEMKSSEP SGSLGLSDYS GLQEELQELR ERYHFLNEEY RALQESNSSL TGQLADLESE RTQRATERWL QSQTLSMTSA ESQTSEMDFL EPDPEMQLLR QQLRDAEEQM HGMKNKCQEL CCELEELQHH RQVSEEEQRR LQRELKCAQN EVLRFQTSHS VTQNEELKSR LCTLQKKYDT SQDEQNELLK MQLQLQTEL QLKVMKSTLV ENQSEKELLC RLQKLHLQH QVTCEKEKLL ERQQQLQEEL QCHEAELQHL RDTVASFKE NEKDTETHAQ LQEMKQLYQA SKDELERQKH MYDQLEQDLL LCQLELKEK ASHPIPEDKG KCANKCDTLL SRLTELQEKY KASQKEMGQL QMEQCELLED QRRMQEEQGQ LQEELHRLTL PLPKSGLLLK SQELLTKLED LCELQLLYQG

image



Product Name: Recombinant Human CCDC136 Protein
Catalog #: PEH5232



Background

Alternative Names

Coiled-coil domain-containing protein 136 (Nasopharyngeal carcinoma-associated gene 6 protein)

Note

For Research Use Only , Not for Diagnostic Use