


Product Name: Recombinant Human RMDN2 Protein
Catalog #: PEH5240



Summary

Name	RMD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96LZ7
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EFLSLGNTFN SITLQDEIHD DQGTTVIFQE RQLQILEKLN ELLTNMEELK EEIRFLKEAI PKLEEYIQDE LGGKITVHKI SPQHRARKRR LPTIQSSATS NSSEEAESG GYITANTDTE EQSFPVPKAF NTRVEELNLD VLLQKVDHLR MSESGKSESF ELLRDHKEKF

image



Background

Alternative Names Regulator of microtubule dynamics protein 2 (RMD-2) (hRMD-2) (Protein FAM82A1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RMDN2 Protein
Catalog #: PEH5240

