

Product Name: Recombinant Human CNTNAP3B Protein
Catalog #: PEH5249



Summary

Name	CNT3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96NU0
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVTFLLSSRSY LALPGNSGED KSVTFQFRT WNRAGHLLFG ELQRGSGSFV LFLKDGKLKL SLFQAGQSPR NVTAGAGLND GQWHSVFSFA KWSHMNVVVD DDTAVQPLVA VLIDSGDTYY FGGCLGNSSG SGCKSPLGGF QGCLRLITIG DKAVDPILVQ QGALGSFRDL

image



Background

Alternative Names Contactin-associated protein-like 3B (Cell recognition molecule Caspr3b)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CNTNAP3B Protein
Catalog #: PEH5249

