


Product Name: Recombinant Human FRG2B Protein
Catalog #: PEH5255



Summary

Name	FRG2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96QU4
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GNCPEKECSL SLNKKRSST PVHNSEIQET CDAHHRGRSR ACTGRSKRHR SRALGVQTPS IRKSLVTSVR AMSEAVYQDL

image



Background

Alternative Names	Protein FRG2-like-1 (FSHD region gene 2 protein family member B) (HSA10-FRG2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use