


**Product Name: Recombinant Human MAGI1 Protein**  
**Catalog #: PEH5256**



## Summary

<b>Name</b>	MAGI1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96QZ7
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	31 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SHSDIVNLIK EAGNTVTLRI IPGDESSNAT LLTNAEKIAT ITTTHTPSQQ GTQETRNTTK PKQESQFEFK APQATQEQDF YTVELERGA K GFGFSLRGGR EYNMDLYVLR LAEDGPAERC GKMRIGDEIL EINGETTKNM KHSRAIELIK NGGRRVRLFL KRGDGSVPEY DPSSDRHGPA TGPQGVPEVR AGPDRRQHPS LESSYPPDLH KSSPHGEKRA HARDPKGSRE YSRQPNEHHT WNGTSRKPDS GACRPKDRAP

## image



## Background

<b>Alternative Names</b>	Membrane-associated guanylate kinase, WW and PDZ domain-containing protein 1 (Atrophin-1-interacting protein 3) (AIP-3) (BAI1-associated protein 1) (BAP-1) (Membrane-associated guanylate kinase inverted 1) (MAGI-1)
--------------------------	--

**Product Name: Recombinant Human MAGI1 Protein**  
**Catalog #: PEH5256**



---

(Trinucleotide repeat-containing gene 19 protein) (WW domain-containing protein 3) (WWP3)

**Note**

For Research Use Only , Not for Diagnostic Use