


Product Name: Recombinant Human CLCC1 Protein
Catalog #: PEH5258



Summary

Name	CLCC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96S66
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSQAKYGISG EKDVPDLSC ADEISECYHK LDSLTYSKIDE CEKKKREDYE SQSNPVFRRY LNKILIEAGK LGLPDENKGD MHYDAEILK RETLLEIQKF LNGEDWKPGA LDDALSDILI

image



Background

Alternative Names	Chloride channel CLIC-like protein 1 (Mid-1-related chloride channel protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use