


Product Name: Recombinant Human TCF12 Protein
Catalog #: PEH5260



Summary

Name	HTF4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99081
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTVVTTEIKT ENKEKDENLH EPPSSDDMKS DDESSQKDIK VSSRGRTSST NEDEDLNPEQ KIEREKERRM ANNARERLRV RDINEAFKEL GRMCQLHLKS EKPQTKLLIL HQAVAVILSL EQQVRERNLN PKAACLKRRE EEKVS AVSAE

image



Background

Alternative Names Transcription factor 12 (TCF-12) (Class B basic helix-loop-helix protein 20) (bHLHb20) (DNA-binding protein HTF4) (E-box-binding protein) (Transcription factor HTF-4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TCF12 Protein
Catalog #: PEH5260

