


Summary

Name	MUC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99102
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IHSTSTSPQE SPAVSQRGHT QAPQTTQESQ TTRSVSPMTD TKTVTTPGSS FTASGHSPSE IVPQDAPTIS AATTFAPAPT GDGHTTQAPT TALQAAPSSH DATLGPSGGT SLSKTGALTL ANSVVSTPGG PEGQWTSASA STSPDTAAAM THTHQAESTE ASGQTQTSEP ASSGSRRTSA

image



Background

Alternative Names	Mucin-4 (MUC-4) (Ascites sialoglycoprotein) (ASGP) (Pancreatic adenocarcinoma mucin) (Testis mucin) (Tracheobronchial mucin) [Cleaved into: Mucin-4 alpha chain (Ascites sialoglycoprotein 1) (ASGP-1); Mucin-4 beta chain (Ascites sialoglycoprotein 2) (ASGP-2)]
--------------------------	--

Product Name: Recombinant Human MUC4 Protein
Catalog #: PEH5261



Note

For Research Use Only , Not for Diagnostic Use