


Summary

| | |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name | K1C12_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q99456 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 38 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | LSGSEKETMQ NLNDRLASYL DKVRALEEAN TELENKIREW YETRGTGTAD ASQSDYSKY PLIEDLRNKI ISASIGNAQL LLQIDNARLA AEDFRMKYEN ELALRQGVEA DINGLRRVLD ELTLTRTDLE MQIESLNEEL AYMKKNHEDE LQSFVRGGPG EVSVEMDAAP GVDLTRLND MRAQYETIAE QNRKDAEAWF IEKSGELRKE ISTNTEQLQS SKSEVTDLRR AFQNLEIELQ SQLAMKKSLE DSLAEAEGDY CAQLSQVQQL ISNLEAQLLQ VRADAERQNV DHQRLLNVKA RLELEIETYR RLLDGEAQGD |

image



Background

Alternative Names Keratin, type I cytoskeletal 12 (Cytokeratin-12) (CK-12) (Keratin-12) (K12)

Product Name: Recombinant Human KRT12 Protein
Catalog #: PEH5263



Note

For Research Use Only , Not for Diagnostic Use