


Product Name: Recombinant Human MNT Protein
Catalog #: PEH5264



Summary

Name	MNT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99583
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LAPAEVKSS EQKKRPGGIG TREVHNKLEK NRR AHLKECF ETLKRNI PNV DDK KTSNLSV LRTALRYIQS LKRKEKEYEH EMERLAREKI ATQQR LAELK HEL SQWMDVL EIDRVLRQTG QPEDDQASTS TASEGEDNID EDM EEDRAGL GPPKLSHRPQ PELLKSTLPP PSTTPAPLPP HPHPHPSVA LPPAHLPVQQ

image



Background

Alternative Names	Max-binding protein MNT (Class D basic helix-loop-helix protein 3) (bHLHd3) (Myc antagonist MNT) (Protein ROX)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MNT Protein
Catalog #: PEH5264

