


Product Name: Recombinant Human NPAS2 Protein
Catalog #: PEH5267



Summary

| | |
|------------------------------|--|
| Name | NPAS2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q99743 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 20 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | RASRNKSEKK RRDQFNVLIK ELSSMLPGNT RKMDKTTVLE KVIGFLQKHN EVSAQTEICD IQQDWKPSFL SNEEFTQLML EALDGFIIAV TTDGSIYVS DSITPLLGHLL PSDVMDQNLL NFLPEQEHSE VYKILSSHML VTDSPSPEYL |

image



Background

Alternative Names Neuronal PAS domain-containing protein 2 (Neuronal PAS2) (Basic-helix-loop-helix-PAS protein MOP4) (Class E basic helix-loop-helix protein 9) (bHLHe9) (Member of PAS protein 4) (PAS domain-containing protein 4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NPAS2 Protein
Catalog #: PEH5267

