


Product Name: Recombinant Human SDF4 Protein
Catalog #: PEH5271



Summary

Name	CAB45_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BRK5
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLNGVKLEMD GHLNRGFHQE VFLGKDLGGF DEDAEPRRSR RKLMMVIFSKV DVNTDRKISA KEMQRWIMEK TAEHFQEAME ESKTHFRAVD PDGDGHVSWD EYKVKFLASK GHSEKEVADA IRLNEELKVD EETQEVLENL KDRWYQADSP PADLLLTEEE FLSFLHPEHS RGMLRFMVKE IVRDLQDGD KQLSVPEFIS LPVGTVENQQ

image



Background

Alternative Names 45 kDa calcium-binding protein (Cab45) (Stromal cell-derived factor 4) (SDF-4)

Note

Product Name: Recombinant Human SDF4 Protein
Catalog #: PEH5271



For Research Use Only , Not for Diagnostic Use