



Summary

Name	TIPIN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BVW5
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLPDYEHVED ETFPPFPPPA SPERQDGEGT EPDEESGNGA PVRVPPKRTV KRNIPKLDAPQ RLISERGLPA LRHVFDKAKF KGKGHEAEDL KMLIRHMEHW AHLRFPKLQF EDFIDRVEYL GSKKEVQTCL KRIRLDLPIL HEDFVSNNDE VAENNEHDVT STELDPFLT N LSEEMFASE LSRSLTEEQQ QRIERNKQLA LERRQAKLLS NSQTLGNDML MNTPRAHTVE EVNTDEDQKE ESNGLNEDIL DNPCNDAIAN TLNEEETLLD QSFKNVQQQL DATSRNITEA

image



Background

Alternative Names TIMELESS-interacting protein

Product Name: Recombinant Human TIPIN Protein
Catalog #: PEH5274



Note

For Research Use Only , Not for Diagnostic Use