


Product Name: Recombinant Human ASCC2 Protein
Catalog #: PEH5288



Summary

Name	ASCC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H1I8
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YYEDEYDDTY DGNQVGANDA DSDELISRR PFTIPQVLRT KVPREGQEED DDDEEDDADE EAPKPDHFVQ DPAVLREKAE ARRMAFLAKK GYRHDSSTAV AGSPRGHGQS RETTQERRKK EANKATRANH NRRTMADRKR

image



Background

Alternative Names Activating signal cointegrator 1 complex subunit 2 (ASC-1 complex subunit p100) (Trip4 complex subunit p100)

Note

For Research Use Only , Not for Diagnostic Use