


Product Name: Recombinant Human UPF3A Protein
Catalog #: PEH5289



Summary

Name	REN3A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H1J1
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEEKRRRREE ERCKKKETDK QKKIAEKEVR IKLLKKPEKG EEPTEKPKE RGEEIDTGGG KQESCAPGAV VKARPMEGSL EEPQETSHSG SDKEHRDVER SQEQESEAQR YHVDDGRRHR

image



Background

Alternative Names Regulator of nonsense transcripts 3A (Nonsense mRNA reducing factor 3A)
(Up-frameshift suppressor 3 homolog A) (hUpf3)

Note

For Research Use Only , Not for Diagnostic Use