


Product Name: Recombinant Human PABPC3 Protein
Catalog #: PEH5295



Summary

Name	PABP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H361
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IFVKNLDKSI N NKALYDTVS AFGNILSCNV VCDENGSKGY GFVHFETHEA AERIAIKKMNG MLLNGRKVFV GQFKSRKERE AELGARAKEF PNVYIKNFGE DMDDERLKDL FGKFGPALSV KVMTDESGKS KGFGFVSFER HEDAQKAVDE MNGKELNGKQ IYVGRAQKKV

image



Background

Alternative Names Polyadenylate-binding protein 3 (PABP-3) (Poly(A)-binding protein 3) (Testis-specific poly(A)-binding protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PABPC3 Protein
Catalog #: PEH5295

