


Product Name: Recombinant Human HEATR1 Protein
Catalog #: PEH5300



Summary

Name	HEAT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H583
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVLARLGDDN IDVVLSAISA FEIFKEHFSS EVTISNLLNL FQRAELSKNG EWYEVLKIAA DILIKEEILS ENDQLSNQVV VCLLPFMVIN

image



Background

Alternative Names	HEAT repeat-containing protein 1 (Protein BAP28) (U3 small nucleolar RNA-associated protein 10 homolog) [Cleaved into: HEAT repeat-containing protein 1, N-terminally processed]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use