


Product Name: Recombinant Human SYNDIG1 Protein
Catalog #: PEH5303



Summary

Name	SYNG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H7V2
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LVHISKISDAG KRNGLINTRN LMAESRDGLV SVYPAPQYQS HRVGASTVPA SLDSSRSEPM QQLDPNTLQ QSVESRYRPN IILYSEGVLR SWGDGVAADC CETTFIEDRS PTKDSLEYPD GKFIDLSADD IKIHTLSYDV EEEEFQELE SDYSSDTESE DNFLMMPPRD

image



Background

Alternative Names	Synapse differentiation-inducing gene protein 1 (SynDIG1) (Dispanin subfamily C member 2) (DSPC2) (Transmembrane protein 90B)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SYNDIG1 Protein
Catalog #: PEH5303

