


Product Name: Recombinant Human ZDHHC11 Protein
Catalog #: PEH5305



Summary

Name	ZDH11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H8X9
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YLINNRKEES SKHQAVRKDP YVQMDKGV LQ QGAGALGSSA QGVKAKSSLL IHKHLCHFCT SVNQDGDSTA REGDEDPCPS ALGAKARNSR LICRRLCQFS TRVHPDGGSM AQEADDAPSI STLGLQQETT EPMKTDSAES

image



Background

Alternative Names Palmitoyltransferase ZDHHC11 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 11) (DHHC-11) (Zinc finger protein 399)

Note

For Research Use Only , Not for Diagnostic Use