


Product Name: Recombinant Human ZMAT3 Protein
Catalog #: PEH5306



Summary

Name	ZMAT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HA38
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAQQAQAHYQ GKNHGGKLRN YYAANSCPPP ARMSNVVEPA ATPVVPVPPQ MGSFKPGGRV ILATENDYCK LCDASFSSPA VAQAHYQGKN HAKRLRLAEA QSNSESESE LGQRRARKEG NEFKMMPNRR NMYTVQNNSA

image



Background

Alternative Names	Zinc finger matrin-type protein 3 (Zinc finger protein WIG-1) (p53-activated gene 608 protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use