



Product Name: Recombinant Human MEAF6 Protein
Catalog #: PEH5307



Summary

Name	EAF6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HAF1
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAMHNKAAPP QIPDTRRELA ELVKRKQELA ETLANLERQI YAFEGSYLED TQMYGNIIRG WDRYLTNQKN SNSKNDRRNR KFKEAERLFS KSSVTSAAAV SALAGVQDQL IEKREPGSGT ESDTSPDFHN QENEPSQEDP EDLDGVSQGV KPQKAASSTS SGGSHSSHKK RKNKNRHRID LKLNKKPRAD Y

image



Background

Alternative Names	Chromatin modification-related protein MEAF6 (MYST/Esa1-associated factor 6) (Esa1-associated factor 6 homolog) (Protein EAF6 homolog) (hEAF6) (Sarcoma antigen NY-SAR-91)
--------------------------	--

Note

Product Name: Recombinant Human MEAF6 Protein
Catalog #: PEH5307



For Research Use Only , Not for Diagnostic Use