


Product Name: Recombinant Human GPBP1L1 Protein
Catalog #: PEH5312



Summary

Name	GPBL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HC44
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LTKLRRRTTD RKSEFLKTLK DDRNGDFSEN RDCDKLEDLE DNSTPEPKEN GEEGCHQNGL ALPVVEEGEV LSHSLEAEHR LLKAMGWQEY PENDENCLPL TEDELKEFHM KTEQLRRNGF

image



Background

Alternative Names Vasculin-like protein 1 (GC-rich promoter-binding protein 1-like 1)

Note

For Research Use Only , Not for Diagnostic Use