


Product Name: Recombinant Human CBX8 Protein
Catalog #: PEH5313



Summary

Name	CBX8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HC52
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERERERERER GTSRVDDKPS SPGDSSKCRG PKPRKELPDP SQRPLGEP SA GLGEYLKGRK LDDTPSGAGK FPAGHSVIQL ARRQDSDLVQ CGVTSPSSAE ATGKLAVDTF PARVIKHRAA

image



Background

Alternative Names Chromobox protein homolog 8 (Polycomb 3 homolog) (Pc3) (hPc3) (Rectachrome 1)

Note

For Research Use Only , Not for Diagnostic Use