


Product Name: Recombinant Human ZNF532 Protein
Catalog #: PEH5314



Summary

Name	ZN532_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCE3
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TLKSIEGPPN LGINLPLSIK PATQNSANQN KEDTKSMNGK EKLEKKSPSP VKKSMETKKV ASPGWTCWEC DCLFMQRDVY ISHVRKEHGK QMKKHPCRQC DKSFSSSHSL CRHNRIKHKG IRKVYACSHC PDSRRTFTKR LMLEKHVQLM HGIKDPLKE MTDATNEEET EIKEDTKVPS PKRKLEEPVL

image



Background

Alternative Names Zinc finger protein 532

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF532 Protein
Catalog #: PEH5314

