


Product Name: Recombinant Human METTL14 Protein
Catalog #: PEH5315



Summary

Name	MET14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCE5
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQKLRRQLLA QQLGAESADS IGAVLNSKDE QREIAETRET CRASYDTSAP NAKRKYLDEG ETDEDKMEYY KDELEMQQDE ENLPYEEIY KDSSTFLKGT QSLNPHNDYC QHFVDTGHRP QNFIRDVGLA DRFEEYPKLR ELIRLKDELI

image



Background

Alternative Names	N6-adenosine-methyltransferase non-catalytic subunit (Methyltransferase-like protein 14) (hMETTL14)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use