


Product Name: Recombinant Human CHD8 Protein
Catalog #: PEH5318



Summary

Name	CHD8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCK8
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGFIWDLISP AENGKTKELQ N HSGLSIPVP RGRKGKKVKS QSTFDIHKAD WIRKYNPDTL FQDESYKKHL KHQCNKVLLR VRMLYYLRQE VIGDQAEKVL GGAIASEIDI WFPVVDQLEV PTTWWDSEAD KLLIGVFKH GYEKYNTMRA DPALCFLEKA GRPDDKAIAA EHRVLDNFSD IVEGVDFDKD CEDPEYKPLQ GPPKDQDDEG DPLMMMDEEI SVIDGDEAQV TQQPGHLFWP

image



Background

Alternative Names Chromodomain-helicase-DNA-binding protein 8 (CHD-8) (EC 3.6.4.12) (ATP-dependent helicase CHD8) (Helicase with SNF2 domain 1)

Note

Product Name: Recombinant Human CHD8 Protein
Catalog #: PEH5318



For Research Use Only , Not for Diagnostic Use