


Product Name: Recombinant Human TUFT1 Protein
Catalog #: PEH5322



Summary

Name	TUFT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NNX1
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HIAQKAGRKT YAMVSSHSAG HSLASELVES HDGHEEIKV YLKGRSGDKM IHEKNINQLK SEVQYIQEAR NCLQKLREDI SSKLDRNLGD SLHRQEIQVV LEKPNGFSQS PTALYSSPPE VDTICINEDVE SLRKTVDLL AKLQEAQRQH QSDCVAFEVT LSRVQREAEQ SNVALQREED RVEQKEAEVG ELQRRLLGME TEHQALLAKV REGEVALEEL RSNNADCQAE REKAATLEKE VAGLREKIH LDDMLKSQQR KVRQMIEQLQ NSKAVIQSKD ATIQLKEKI AYLEAENLEM

image



Background

Alternative Names Tuftelin

Product Name: Recombinant Human TUFT1 Protein
Catalog #: PEH5322



Note

For Research Use Only , Not for Diagnostic Use