


Product Name: Recombinant Human MED4 Protein
Catalog #: PEH5328



Summary

Name	MED4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NPJ6
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KERLGGGLGV AGGNSTRERL LSALEDLEVL SRELIEMLA SRNQKLLQAG EENQVLELLI HRDGEFQELM KLALNQGKIH HEMQVLEKEV EKRDSIDIQQL QKQLKEAEQI LATAVYQAKE KLKSIEKARK

image



Background

Alternative Names	Mediator of RNA polymerase II transcription subunit 4 (Activator-recruited cofactor 36 kDa component) (ARC36) (Mediator complex subunit 4) (TRAP/SMCC/PC2 subunit p36 subunit) (Vitamin D3 receptor-interacting protein complex 36 kDa component) (DRIP36)
--------------------------	--

Note

Product Name: Recombinant Human MED4 Protein
Catalog #: PEH5328



For Research Use Only , Not for Diagnostic Use