


Product Name: Recombinant Human EQTN Protein
Catalog #: PEH5329



Summary

Name	EQTN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQ60
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KPTIEALPNV LPLNEDVKNQ EEKNEDHTPN YAPANEKNGN YYKDIKQYVF TTQNPNGTES EISVRATDDL NFALKNDKTV NATTYEKSTI EEETTSEPS HKNIQRSTPN VPAFWTMLAK AINGTAVVMD DKDQLFHPIP ESDVNATQGE NQPDLLEDLKI

image



Background

Alternative Names Equatorin (Acrosome formation-associated factor)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EQTN Protein
Catalog #: PEH5329

