


Summary

| | |
|------------------------------|--|
| Name | A1CF_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9NQ94 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 65 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | LSGTQKEAAL RALVQRTGYS LVQENGQRKY GGPPPGWDAA PPERGCEIFI GKLPRDLFED ELIPLCEKIG KIYEMRMMMD FNGNNRGYAF VTFSNKVEAK NAIKQLNNYE IRNGRLLGVC ASVDNCRFLV GGIPKTKKRE EILSEMKKVT EGVVDVIVYP SAADKTKNRG FAFVEYESHR AAAMARRKLL PGRIQLWGHG IAVDWAEPEV EVDEDTMSSV KILYVRNLML STSEEMIEKE FNNIKPGAVE RVKKIRDYAF VHFSNREDAV EAMKALNGKV LDGSPIEVTL AKPVDKDSYV RYTRGTGGRG TMLQGEYTYS LGQVYDPTTT YLGAPVFYAP QTYAAIPSLH FPATKGHLSN RAIIRAPSVR EIYMNVPVGA AGVRLGGRG YLAYTGLGRG YQVKGDKRED KLYDILPGME LTPMNPVTLK PQGIKLAPQI LEEICQKNNW GQPVYQLHSA IGQDQRQLFL YKITIPALAS QNPAIHPFTP PKLSAFVDEA KTYAAEYTLQ TLGIPTDGGD GTMATAAAAA TAFPGYAVPN ATAPVSAAQL KQAVTLGQDL AAYTTYEVYP |

image



Product Name: Recombinant Human A1CF Protein
Catalog #: PEH5330



Background

Alternative Names APOBEC1 complementation factor (APOBEC1-stimulating protein)

Note

For Research Use Only , Not for Diagnostic Use