


Product Name: Recombinant Human RPRD1B Protein
Catalog #: PEH5333



Summary

Name	RPR1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQG5
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEKKSLKRTF QQIQEEEDDD YPGSYSPQDP SAGPLLTEEL IKALQDLENA ASGDATVRQK IASLPQEVQD VSLLEKITDK EAAERLSKTV DEACLLAEY NGRLAAELED RRQLARMLVE YTQNQKDVLS EKEKKLEEKY QKLARVTQVR KELKSHIQSL PDLSELLPNVT

image



Background

Alternative Names Regulation of nuclear pre-mRNA domain-containing protein 1B (Cell cycle-related and expression-elevated protein in tumor)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RPRD1B Protein
Catalog #: PEH5333

