


Product Name: Recombinant Human DROSHA Protein
Catalog #: PEH5338



Summary

Name	RNC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRR4
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELGSRQEKAK AARPPWEPPK TKLDEDLESS SESECESDSE STCSSSDSE VFDVIAEIKR KKAHPDRLHD ELWYNDPGQM NDGPLCKCSA KARRTGIRHS IYPGEEAIKP CRPMTNNAGR

image



Background

Alternative Names	Ribonuclease 3 (EC 3.1.26.3) (Protein Drosha) (Ribonuclease III) (RNase III) (p241)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use