


Product Name: Recombinant Human KRT82 Protein
Catalog #: PEH5340



Summary

Name	KRT82_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NSB4
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEQIKCLNNR FASFINKVRF LEQKNKLLT KWNFMQQQRC CQTNIEPIFE GYISALRRQL DCVSGDRVRL ESELCSLQAA LEGYKKKYEE ELSLRPCVEN EFVALKKDVD TAFLMKADLE TNAEALVQEI DFLKSLYEEE ICLLSQISE TSVIVKMDNS RELDVDGIIA EIKAQYDDIA SRSKAEAEAW YQCRYEELRV TAGNHCDNLR NRKNEILEMN KLIQRLQQET ENVKAQRCKL EGAIAEAEQQ GEAALNDAKC KLAGLEEALQ KAKQDMACLL KEYQEVMSK LGLDIEIATY RRLLEGEEHR

image



Background

Alternative Names	Keratin, type II cuticular Hb2 (Keratin-82) (K82) (Type II hair keratin Hb2) (Type-II keratin Kb22)
--------------------------	---

Product Name: Recombinant Human KRT82 Protein
Catalog #: PEH5340



Note

For Research Use Only , Not for Diagnostic Use