


Product Name: Recombinant Human PDS5B Protein
Catalog #: PEH5345



Summary

Name	PDS5B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NTI5
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGRLDSSEMD HSENYDTMS SPLPGKSDK RDDSDLVRSE LEKPRGRKKT PVTEQEEKLG MDDLTKLVQE QKPKGSQRSR KRHTASESD

image



Background

Alternative Names	Sister chromatid cohesion protein PDS5 homolog B (Androgen-induced proliferation inhibitor) (Androgen-induced prostate proliferative shutoff-associated protein AS3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use