


Product Name: Recombinant Human INKA2 Protein
Catalog #: PEH5346



Summary

Name	INKA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NTI7
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VCGRDLAPLP RTQPHQSCAQ QGPERVEPDD WTSTLMSRGR NRQPLVLGDN VFADLVGNWL DLPELEKGGE KGETGGAREP

image



Background

Alternative Names	PAK4-inhibitor INKA2 (Induced in neural crest by AP2-alpha protein-related homolog) (HInca-r) (Inka box actin regulator 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use